



Introduction to QA/QC in analytical radiochemistry

Marko Štrok Jožef Stefan Institute



This project has received funding from the European Union's Horizon 2020 research and innovation programme under the Euratom research and training programme 2014-2018 under grant agreement No 754 972





QA/QC

• Crucial for early detection of any problems with instrument orm ethod

- Enable that we can take appropriate measures to correct for them
- Accreditation to ISO 17025





Method validation

• Fit for purpose

- Analysis of reference materials
- Participation in interlaboratory comparisons
- Linearity
- Reproducibility
- Repeatability
- Validated methods need to be checked before used in our laboratory






Parameters affecting quality of results

- Room temperature
- Detection efficiency
- Instrument background
- Blank count rate
- Results of analysis of reference materials or internal standards
- Plot them with time => control chart
- At least 10 inputs, better through whole year





Control chart








Using control chart

- Method is under control
 - value within upper and lower warning limits or
 - between warning and action limit and last two values are within upper and lower warning llimits
- Method is under control but out of statistical control
 - 7 consecutive values are increasing or decreasing
 - 10 out of eleven last consecutive values lies on one side of average value but within upper and lower limits
- Method is not under control:
 - control value is out of action limit
 - control value is betweeen warning and action limit and at least one out of last two values are is also between warning and action limit







Institut "Jožef Stefan", Ljubljana, Slovenija



This project has received funding from the European Union's Horizon 2020 research and innovation programme under the Euratom research and training programme 2014-2018 under grant agreement No 754 972