



Basics of analytical radiochemistry

Marko Štrok Jožef Stefan Institute




This project has received funding from the European Union's Horizon 2020 research and innovation programme under the Euratom research and training programme 2014-2018 under grant agreement No 754 972






Analytical radiochemistry and radioactivity

- Analytical chemistry dealing with radioactive elements or radionuclides
- Radionuclide = nuclide, which has excess nuclear energy making it unstable
- Alpha, beta, gamma emitters






- We can use radioacivity to detect them
- Radionuclides are present everywhere






Radioactive decay

• Radionuclide can decay and/or gow up



• Radioactive equilibrium








Natural and anthropogenic radionuclides


• Anthropogenic = produced by human activities










Natural = produced by other activities













Natural radionuclides

Primordial = half life long enough to not decay since creation of Earth ulletU-238, U-235, Th-232, K-40







- Secondary = decay products of primordial, but with shorter half life Ra-226, Pb-210, Po-210,...
- Cosmogenic = produced with interaction with cosmic rays H-3, C-14, Be-7, Be-10...







Typical analysis steps in analytical radiochemistry

- **Sample pre-treatment** (filtration, sieving, homogenisation, digestion, dissolution, pre-concentration,...)
- **Separation of analyte** (precipitation, liqud extraction, extraction chromatography, spontaneous deposition, distillation,...)
- Preparation of counting source (micro-coprecipitation, electrodeposition, precipitation, mixing with scintillation cocktail,...)
- **Measurement** (radiometry, mass spectrometry,...)
- Data treatment and calculation of results





Concept of tracers in analytical radiochemistry

- To trace lossess during analytical procedure
- Ideal tracer should:
 - behave the same as our analyte
 - not be present in our sample
 - have sufficiently long half live
 - have known concentration traceable to SI units
 - be readily available
- Some examples of tracers:

U-232 (uranium), Th-229 (thorium), Ba-133 (radium), Po-209 (Po-210), stable Pb (Pb-210), stable Sr or Sr-85 (Sr-90)





Radiometry


Gamma spectrometry

Alpha spectrometry



Liquid scintillation counting Proportional counting










Mass spectrometry









Institut "Jožef Stefan", Ljubljana, Slovenija



This project has received funding from the European Union's Horizon 2020 research and innovation programme under the Euratom research and training programme 2014-2018 under grant agreement No 754 972